

RAD Shoulder System Instructions/ Instrucciones por el RAD Sistema del Hombro

To Put On the RAD Shoulder System Brace/Para poner el RAD Sistema del Hombro:

1. Position your arm in the RAD Shoulder system arm sling./ *Posición el brazo en el bolsillo del brazo del RAD Sistema del Hombro.*

2. Secure the waist strap using the buckle attached at the front of the lateral waist support. /*Abrochar la correa a la hebilla*



3. Put the strap with the RAD label on it around your shoulder as if you were putting on a back pack. You should be able to see the letters RAD on your shoulder when you look in the mirror./ *Ponga la correa con las letras RAD por la espalda y alrededor del hombro como ponerse una mochila. Usted debe poder de leer las letras RAD en un espejo.*



4. Secure the loose strap attached to the RAD Shoulder System three point harness to using the hook to attach at the forearm of the arm sling./ *Asegurar la correa al gancho en la parte delantera y cerca de la mano.*



5. Attach the strap on the sling over your arm. / *Asegurar la correa con el Velcro sobre el brazo.*



To Take Off the RAD Shoulder System Brace/Para quitar el RAD Sistema del Hombre:

1. Detach the strap located over the arm sling./ *Sacar la correa situada sobre el brazo.*



2. Detach the strap at the forearm of the sling./*Sacar la correa cerca de la mano al gancho.*



3. Unbuckle the waist strap. /*Desabrochar la correa a la hebilla.*



4. Take the strap with the RAD label on it off of your shoulder and slide the RAD Shoulder system off of your arm. /*Quitar la correa con las letras RAD por la espalda y quitar el RAD Sistema del hombro por el brazo con cuidado.*



How to use the RAD Shoulder System Brace as a sling /*Para usar el RAD Sistema del Hombre con sólo el bolsillo del brazo:*

- 1.** Position your arm in the RAD Shoulder System arm sling. /*Posicione el brazo en el bolsillo del brazo del RAD Sistema del Hombre.*
- 2.** Pass the RAD Shoulder System three point harness behind the back and over the uninjured shoulder by placing your free arm through the ring created by the straps. / *Ponga la correa por la espalda y alrededor del hombro.*
- 3.** Secure the loose strap attached to the RAD Shoulder System three point harness to the d-ring at the forearm of the arm sling. / *Cierra la correa de la cadena al frente.*

